

by using the UR-RTDE library developed by SDU [88]. It provides a real-time interface to exchange I/O data and control the robot arm from an external application or program. The library works by establishing a communication link between the computer and the UR robot arm using an Ethernet connection. It utilizes the Real-Time Data Exchange (RTDE) protocol, which is a proprietary protocol developed by Universal Robots for real-time communication [89]. The connection makes it possible to receive real-time feedback from the robot arm, such as joint positions, joint velocities, and other sensor data. Additionally, it enables sending commands and controlling the robot arm's motion in real-time.

The pneumatic clamps and the cylinder are controlled by using electrically actuated 5/2-way control valves. These are both connected to two digital output connectors of the UR control box. Each output can control the airflow in one output 6 mm tube of the valve, which is responsible for extending or compressing the clamps and the cylinder. The Robotiq Hand-E gripper is connected to the robot arm by a power and communication cable with a USB adapter. The microscope, camera, and rotary electric gripper are all connected by a USB cable (type 3.0 output) to the PC. While the camera and the gripper can be controlled by using specific Python libraries (pyrealsense2, minimalmodbus), the microscope needs the CytoSMART driver installed on the PC. The entire system is controlled by an application written in Python. This includes a simple human-machine interface to control the tasks and see the output of the microscope. The recorded data (images taken by microscope and camera) are stored locally on the PC. The functional interconnection of hard- and software is visualized in fig. 3.6. The UR control box, the rotary electric gripper, the incubator, and the PC are connected to standard 230V electrical sockets.

3.1.3 Process and workflow description

Similar to the manual process, the autonomous process is split into three high-level workflows.

- Workflow A: Analyzing cell growth.
- Workflow B: Changing media.
- Workflow C: Passaging.

The system should be able to execute the three workflows independently, depending on the user input or the scheduled plan. To enhance reusability and simplicity, the steps to achieve completion of the workflows are structured into workflow modules.

The following gives a mid-level overview of the purpose of each process module, including a summarized description of the substeps. A low-level description is given in THE APPENDIX???

- Module 1: Get a cell culture flask from the incubator (open clamps, open door, grip the flask, move the flask outside, close door and clamps)
- Module 2: Analyze cells (move flask to regripping station, regrip, place the flask on the microscope, take images at different positions, move back to start position)
- Module 3: Place a flask in a flask holder (move to the flask holder, place the flask inside, and move back to the start position)
- Module 4: Decap flask(s) (move to the flask, grip the lid, take off the lid, place the lid on a lid holder, and move back to the start position)
- Module 5: Remove liquid from a flask (take the open flask from a flask holder, move to the waste container, pour out all the liquid, and place the flask back into the flask holder)
- Module 6: Add liquid to a flask (heat up media or washing solution, take the bottle, decap the bottle, regrip the bottle, pour a specific amount into the flask, regrip the bottle, cap the




Figure 3.6: Simplified overview of the functional interconnection of hard- and software.

bottle, and place it back into the bottle holder)

- Module 7: Add trypsin to a flask (move to the trypsin unit, move up the bottle dispenser, push it down, and move back to the start position)
- Module 8: Get three empty flasks and place them into the flask holders (move to flask storage, grip an empty flask, move to the flask holder, place the flask inside, and move back to the start position (3 times))
- Module 9: Split cells into empty flasks (take the full flask from the flask holder, move to the empty flasks, pour liquid three times, and place the empty flask into the flask storage)
- Module 10: Cap flask(s) (get a lid from the lid holder, place the lid on a flask, cap the flask, and move back to the start position)
- Module 11: Place a flask in the incubator (open clamps, open door, move the flask inside, place the flask on the flask storage, close door, and close clamps)

